




Speciality Uses for the COLLIS CURVE TOOTHBRUSH



Dental Braces

Collis Curve Medium and Periodontal brushes are recommended for braces. The curved bristles work along under the ligatures connecting the arch wires to remove materia alba. Discolouration due to calcification along the gum line can be minimised by proper cleaning. With a minimum of placement technique the curved bristle will work its way around the bands and remove debris along the gum line. Brush horizontally allowing curved bristles to enter between bands and wires to remove debris caught in brackets.



Dental Implants

The Collis Curve Periodontal brush is also known as our implant brush. It is unique in its ability to reach up even higher into the sulcular areas associated with implants. Since with implants there is no periodontal membrane surrounding the implant/abutment, there is very little protection against the invasion of virulent bacteria. The patient is really dealing with pockets rather than sulci. Antibacterial mouth washes can be introduced into sulcular spaces by simply wetting the brush and letting the curved bristles enter the sulcus and carry along the agent. With the Collis Curve toothbrush there is a return to **simplicity** in maintenance.

Simplicity in **Implant/Abutment** supported fixed prostheses. With a minimum of placement technique the curved bristles will work their way into the artificial sulcus created by the emergency profile of the implant to remove debris

Simplicity in **Bar/Implant** supported full or partial dentures. The natural curve of the Collis Curve implant brush will surround and clean the bar and cantilever introducing antibacterial solutions via the wet, curved bristles that work into sulcular spaces.

Telephone: 01934 862415
Email: CollisCurveUK@aol.com
www.colliscurve.co.uk

