

Toothbrush Selection

The five different styles are denoted by the colour of the cap on the tube.



RED TOP - PERIODONTAL


CLEAR TOP - SOFT

BLACK TOP - MEDIUM

BLUE TOP - JUNIOR

GREEN TOP - BABY

PERIO



The Collis Curve Perio brush is for patients with recessed gums, long or larger teeth, patients with dental appliances such as braces, caps, crowns, implants etc. And for patients who have periodontal pockets. The longer, curved bristles reach gently up around these special areas. Bristle diameter is 0.009" with longer trim height; adult length 1/2". See additional material for implant and orthodontic applications.

MEDIUM



The Collis Curve Medium bristle is the 'stiffest' we have. It will give you a positive feeling of massage and do a superior job of cleaning. This brush is good to recommend for brushers who have healthy gums and have already used the Collis Curve previously. Also for cleaning braces. Bristle is 0.009" trim height; adult length 15/32". See additional material on the medium brush for orthodontic appliances.

SOFT



The Collis Curve soft brush is a gentle brush designed to increase circulation to the gum line area by its lightly massaging action. Ideal for carers and brushers with a limited range of motion. Indications for use include: sensitive gums; gums bleed upon brushing; post-operative conditions, and a first time user of the curved bristle. This is a good brush to use until the gums are in a healthy condition at which point a medium bristle will give you a more positive feeling of massage. Bristle diameter 0.007" trim height,; adult length 15/32"

JUNIOR



The Collis Curve Junior brush has a small head for easy manipulation in the confined spaces of the young person's mouth. Start using at the appearance of the first molar, usually at 6 years old. It is recommended to assist children with brushing until at least 6 years of age.

BABY



The Collis Curve Baby has a smaller head for easy manipulation by parents as they brush their children's teeth. Collis Curve's unique 'on track' comfortable back and forth technique makes it easy for parents to brush their children's teeth. Soft bristles 0/007", small brush head. Use at appearance of first tooth

www.colliscurve.co.uk

Telephone: +44(0)1934 862415

Email: colliscurveuk@aol.com

1 Withey Cottages, Blackmoor, Langford, BS40 5HW, United Kingdom